

molbox1 MFC Interface Connector (25 pin DSUB) Pin No.	Signal Description
1	Case Ground
2	Supply Common
3	MFC Output (+)
4	+15 Volts
5	N/C
6	N/C
7	Set Point Sense (-)
8	N/C
9	N/C
10	N/C
11	Valve Test (-)
12	Current MFC output (+)
13	Current MFC Set (+)
14	MFC Setpoint (+)
15	MFC Setpoint (-)
16	MFC Output (-)
17	Valve Test (+)
18	N/C
19	-15 Volts
20	Set Point Sense (+)
21	N/C
22	N/C
23	N/C
24	RS485 TXRX+
25	RS485 TXRX-